



*Bringing Reliability, Innovation
and Value to Health Care*

Elixir Surgical is much more than a supplier of high precision surgical instruments: a partner who listens carefully and contributes ideas, who understands the customers' needs and requirements and provides for them on a high level so that they can perform at the highest level.



Elixir Surgical manufactures a complete range of gynaecology instruments including but not limited to vaginal speculums, endocervical speculums, uterine dilators, uterine sounds, uterine depressors, uterine polypus forceps, uterine seizing forceps, uterine curettes, tenaculum forceps, myoma instruments, fistula instruments, cervical biopsy punches and IUD instruments.

Elixir Surgical

Web: www.elixirsurgical.com E-mail: info@elixirsurgical.com

OB/GYN INSTRUMENTS



770212

OB Delivery Forceps

ELLIOTT OBSTETRICAL FORCEPS

Set screw in handle
Weck-Pattern, Bright Finish

770212 15 1/2" (39.3 cm)



214560

ELLIOTT OBSTETRICAL FORCEPS

Pilling-Pattern, Satin Finish

214560 13" (33 cm)



214575

KIELAND OBSTETRICAL FORCEPS

214575 14 1/4" (36.1 cm)



770233

LUIKHART OBSTETRICAL FORCEPS

770233 15 1/2" (39.3 cm)



214600

McLEAN-TUCKER OBSTETRICAL FORCEPS

214600 15 1/2" (39.4 cm)



214605

McLEAN-TUCKER-LUIKART OBSTETRICAL FORCEPS

214605 15 1/2" (39.4 cm)

MILLIMETERS

OB Delivery Forceps (continued)



770320

PIPER OBSTETRICAL FORCEPS

770320 18" (45.7 cm)



214550

SIMPSON OBSTETRICAL FORCEPS

214550 14" (35.5 cm)

214555 12" (30.4 cm)



214570

SIMPSON-LUIKART OBSTETRICAL FORCEPS

214570 14" (35.5 cm)



214615

BILL TRACTION HANDLE

For Obstetrical Forceps

214615 7 1/2" (19 cm)



214595

DE LEE OBSTETRICAL FORCEPS

214595 12" (30.4 cm)

214596 14 1/4" (36.1 cm)



214635

DEWEY, AXIS TRACTION

214635 17 3/4" (45 cm)



INCHES

OB/GYN INSTRUMENTS



770237

OB Delivery Forceps *(continued)*

MURLESS FETAL HEAD EXTRACTOR

770237 10" (25.4 cm)



214770

OB Miscellaneous Delivery

BETHEA SHEET HOLDERS

214770 1" (2.5 cm) wide, 5 3/8" (12.9 cm) diameter



834360

DE LEE HILLIS FETAL STETHOSCOPE

834360



214090

UMBILICAL SCISSORS

214090 4 1/8" (10.4 cm)

OB Miscellaneous Delivery *(continued)*



770729

HASSELLTINE UMBILICAL CLAMP

2 3/8" (6 cm)

Package of 12

770729



742110

GOMCO CIRCUMCISION CLAMP

742100	Extra Small	11 mm
742102	Newborn	13 mm
742104	Infant	5 mm
742106	Child	16 mm
742108	Youth	21 mm
742110	Adult	26 mm



212001

GYN Retractors, Specula

GRAVES VAGINAL SPECULA

212000	Small. 3 9/32" x 7/8" blades (8.3 cm x 2.2 cm)
212001	Medium. 4 1/4" x 1 3/8" blades (10.7 cm x 3.4 cm)
212002	Large. 5" x 1 15/32" blades (12.7 cm x 3.7 cm)
212005	Right open-sided Vaginal Specula, Medium. 4 3/8" x 1 1/32" (11.1 cm x .8 cm)



212005



INCHES

OB/GYN INSTRUMENTS



212011

GYN Retractors, Specula *(continued)*

PEDERSON VAGINAL SPECULA

Narrow, straight side blades

- 212010 Small. 3" 3" x 1/2" (7.6 cm x 1.3 cm) blades
- 212011 Medium. 4 3/8" x 3 1/32" (11.1 cm x 2.4 cm) blades
- 212012 Large. 4 3/4" x 1 1/32" (12 cm x 2.6 cm) blades



212151

GUTTMAN OBSTETRICAL RETRACTORS

For post-partum inspection and repair. Swivel blades. Maximum spread is 2" x 4" (5 x 10 cm)

- 212150 Standard 3 7/8" (9.8 cm) blades
- 212151 Extra large 4 3/16" (10.6 cm)



212196

AUVARD WEIGHTED VAGINAL SPECULUM

- 212195 2 3/4" (6.9 cm), weight 2 3/4 lbs. (1.2 kg)
- 212196 3 1/8" (7.9 cm), weight 2 3/4 lbs. (1.1 kg)
- 212197 4" (10 cm), weight 3 lbs. (1.4 kg)



212185

GARRIGUE WEIGHTED VAGINAL SPECULUM

- 212185 Cut-out handle allows drainage. Blade 1 21/32" x 3" (4.2 cm x 7.6 cm) Weight 2 lbs. Chrome plated



212187

GYN Retractors, Specula *(continued)*

HENROTIN WEIGHTED VAGINAL SPECULUM

- 212187 Blade 1 3/8" x 3 1/4" (3.4 cm x 8.3 cm)
Weight 3 lbs. Chrome plated



212227

SIMS VAGINAL RETRACTORS

Double-end

- 212225 7/8" x 2 5/8" and 15/16" x 3 1/8" (2.2 cm x 6.6 cm) and (2.3 cm x 7.9 cm)
- 212226 1" x 2 1/2" and 1 1/4" x 3" (2.5 cm x 6.3 cm) and (4.4 cm x 7.5 cm).
- 212227 15/16" x 3 1/4" x 1 3/16" x 3 1/2" (2.3 cm x 8.2 cm) and (3 cm x 8.8 cm)
- 212228 1 1/4" x 3" and 1 1/2" x 3 1/4" (3.2 cm x 7.5 cm) and (3.8 cm x 8.3 cm).
- 212229 1 3/16" x 3 1/4" x 1 7/16" x 4" (3 cm x 8.2 cm) and (3.6 cm x 10.1 cm)
- 212230 1 7/16" x 3 1/4" x 1 5/8" x 4" (3.6 cm x 8.2 cm) and (4.1 cm x 10.1 cm).



212241

SIMS VAGINAL RETRACTORS

Single-end

- 212241 1 1/8" x 3 5/8" (2.8 cm x 9.2 cm)
- 212242 1 3/8" x 3 1/2" (3.4 cm x 8.8 cm)



212264

JACKSON VAGINAL RETRACTORS

- 212263 1 9/16" x 3" (4 cm x 7.5 cm)
- 212264 1 9/16" x 3 19/32" (3.9 cm x 9.1 cm)
- 212265 1 19/32" x 4" (4 cm x 10.2 cm)



212274

EASTMAN VAGINAL RETRACTORS

- 212271 1 3/8" x 3 1/2" (3.5 cm x 8.8 cm)
- 212274 1 9/16" x 3" (4 cm x 7.5 cm)
- 212275 1 9/16" x 3 1/2" (4 cm x 9 cm)

OB/GYN INSTRUMENTS

GYN Retractors, Specula (continued)

HEANEY HYSTERECTOMY RETRACTORS

- 212300 1 1/8" x 1 3/4" (2.8 cm x 4.4 cm) blade
- 212301 1 1/8" x 2 7/8" (2.8 cm x 7.4 cm) blade
- 212302 1 1/8" x 3 5/16" (2.8 cm x 8.4 cm) blade

HEANEY-SIMON VAGINAL RETRACTOR

- 212305 1" x 4 3/8" (2.5 cm x 11.1 cm)
Blade length 11 1/2" (29.2 cm)

PRYOR-PEAN RETRACTOR

- 212295 10 1/2" (26.6 cm) long
Blade 4 7/16" x 1" (11.2 cm x 2.5 cm)

SHORT MODEL GELPI, SELF-RETAINING RETRACTOR

- 165439 5 1/4" (13.3 cm) overall length
2 1/2" (6.3 cm) maximum spread

FACE LIFT RETRACTOR

- 525215 Fiberoptic lighting, light exits on
top of Deaver type retractor.
11" (28 cm)

212301

212305

212295

165439

525215





165440

GYN Retractors, Specula (continued)

GELPI PERINEAL RETRACTOR

Single sharp prongs. Improved ratchet mechanism. Stainless steel. 7" (17.8 cm) overall length.

165440 Standard type (no stops)



485410

O'SULLIVAN-O'CONNOR VAGINAL RETRACTORS

Spread 52.0 mm

Blades on Frame 52 x 10 mm

Blades Detachable 4 1/4" (10.7 cm)

Small 17 x 45 mm

Medium 37 x 88 mm

Large 43 x 90 mm

485410 6 3/4" (17.1 cm)



212350

DE LEE UNIVERSAL RETRACTOR

212350 2 3/4" x 2" (7 cm x 5 cm) blade



INCHES

OB/GYN INSTRUMENTS



164522

GYN Retractors, Specula (continued)

DOYEN HYSTERECTOMY RETRACTORS

Concave blade. Length 9 1/2" (24.1 cm)

- 164520 1 3/4" x 2 1/2" (4.4 cm x 6.3 cm) blade
- 164521 2 1/8" x 3" (5.4 cm x 7.7 cm) blade
- 164522 2 1/4" x 3 1/2" (5.8 cm x 9.0 cm) blade



884227

GYN Needle Trumpet, Tumor Screw and I.U.D. Hook

IOWA TRUMPET

5 5/8" (14.2 cm)
Needle Guide for Paracervical Anesthesia
Ball Tip 7 mm

884227



752700

DOYEN TUMOR SCREW

Spiral 36 x 25 mm

752700 5 3/4" (14.6 cm)



752529

I.U.D. HOOK

752529 11 1/16" (28 cm)



746799

GYN Catheter Guides

CATHETER GUIDE

15" (38.1 cm)
See Reference Section for French Size/mm Conversion Chart

746799 9 French



740165

FEMALE CATHETER SILVER

- 740163 12 Fr.
- 740165 14 Fr.
- 740166 16 Fr.
- 740167 18 Fr.



GYN Sounds and Dilators



212600

EMMETT UTERINE PROBE

212600 11" (28 cm), Stainless steel



212605

SIMPSON UTERINE SOUND

212605 Graduated in inches.
13 1/4" (33.6 cm), Silver plated



212610

SIMS UTERINE SOUNDS

212610 Non-graduated, 12" (30 cm)
Silver plated
212612 Graduated in cm.
13 1/4" (33.6 cm)



740223

VAN BUREN URETHRAL SOUND

11" (27.9 cm)
Curved 15 degree
See Reference Section for French Size/mm Conversion Chart

740223 34 French



746772

JEWETT SOUNDS

11 1/2" (29.2 cm)
See Reference Section for French Size/mm Conversion Chart

746772 12 French
746773 14 French
746774 16 French



212668

HANK UTERINE DILATORS

Double-ended. 11 sizes. 9-10 Fr. to 29-30 Fr.
Length 10 3/4" (27.3 cm). Single dilator only on sizes
9-10 Fr. to 29-30 Fr.
Single Units
Pilling-Pattern, Satin Finish

212661 9-10 Fr. (3.0 mm x 3.3 mm)
212662 11-12 Fr. (3.6 mm x 4.0 mm)
212663 13-14 Fr. (4.3 mm x 4.7 mm)
212664 15-16 Fr. (5.0 mm x 5.3 mm)
212665 17-18 Fr. (5.7 mm x 6.0 mm)
212666 19-20 Fr. (6.3 mm x 6.7 mm)
212667 21-22 Fr. (7.0 mm x 7.3 mm)
212668 23-24 Fr. (7.7 mm x 8.0 mm)
212669 25-26 Fr. (8.3 mm x 8.7 mm)
212670 27-28 Fr. (9.0 mm x 9.3 mm)
212671 29-30 Fr. (9.3 mm x 10 mm)
212659 Hank's Uterine Dilators, Set of 6
(French Scale) 11-12 Fr. to 21-22 Fr.
(3.66 mm to 7.33 mm)

OB/GYN INSTRUMENTS

GYN Sounds and Dilators *(continued)*

HANK UTERINE DILATORS

10 1/2" (26.6 cm) American Scale
End Lengths 3" (7.6 cm)
Double Ended
Weck-Pattern, Bright Finish



752382

- 752380 Size 9/10
- 752381 Size 11/12
- 752382 Size 13/14
- 752383 Size 15/16
- 752384 Size 17/18
- 752385 Size 19/20

HEGAR UTERINE DILATORS

8" (20.3 cm)
Double Ended, Hollow
Weck-Pattern, Bright Finish



752354

- 752338 Size 1/2 mm
- 752347 Size 19/20 mm
- 752348 Size 21/22 mm
- 752349 Size 23/24 mm
- 752351 Size 25/26 mm
- 752354 Complete Set in Case

HEGAR UTERINE DILATORS

Double-ended, length 8" (20.3 cm)
Pilling-Pattern, Satin Finish

- 212673 Set of 8, 3 mm to 18 mm

HEGAR UTERINE DILATORS

Double-ended
Pilling-Pattern, Satin Finish



212681

- 212675 3 mm to 4 mm
- 212676 5 mm to 6 mm
- 212677 7 mm to 8 mm
- 212678 9 mm to 10 mm
- 212679 11 mm to 12 mm
- 212680 13 mm to 14 mm
- 212681 15 mm to 16 mm
- 212682 17 mm to 18 mm



MILLIMETERS



212647

GYN Sounds and Dilators *(continued)*

PRATT UTERINE DILATORS

Double-ended. 8 sizes. 13-15 Fr. to 41-43 Fr.
Length 11 1/2" (29.2 cm), Stainless steel
Pilling-Pattern, Satin Finish

212647 Set of 8 in canvas case



752394

PRATT UTERINE DILATORS

11 1/2" (29.2 cm)
Double Ended
Weck-Pattern, Bright Finish

- 752390 Size 13/15 French
- 752391 Size 17/19 French
- 752392 Size 21/23 French
- 752393 Size 25/27 French
- 752394 Size 29/31 French
- 752395 Size 33/35 French
- 752396 Size 37/39 French
- 752397 Size 41/43 French



212715

GOODELL UTERINE DILATORS

- 212715 Small, 11 1/4" (28.6 cm)
- 212716 Large, 13" (33 cm)



213120

GYN Tenaculum

ADAIR TENACULUM

213120 7" (17.8 cm)



213110

BARRETT TENACULUM

213110 7" (17.8 cm)



INCHES

OB/GYN INSTRUMENTS



GYN Tenaculum *(continued)*



213100

BRAUN TENACULUM

213100 9 1/2" (24.1 cm)



754106

BRAUN UTERINE TENACULUM

Curved Sideways

754106 9" (22.8 cm)



213090

STAUDE TENACULUM

213090 5 1/2" (14 cm)



213095

STAUDE-MOORE TENACULUM

213095 6 3/8" (16.1 cm)



754160

WYLIE UTERINE TENACULUM

1 x 2 Teeth

754160 9 1/2" (24.1 cm)

GYN Tenaculum *(continued)*



213085

DUPLAY TENACULUM

213085 10" (25.4 cm)



754125

HULKA TENACULUM

10 3/4" (27.3 cm)
Probe 5" (12.7 cm) with 3.5 mm Tip
Single Tooth

754125



213065

JARCHO TENACULUM

213065 7 3/4" (19.6 cm)



213105

KAHN TENACULUM

213105 9 1/2" (24.3 cm)



213080

SCHROEDER TENACULUM

213080 10" (25.4 cm)



INCHES

OB/GYN INSTRUMENTS



213182

GYN Vulsellum

SCHROEDER VULSELLUM

- 213180 Straight, 9 1/2" (24.1 cm)
- 213182 Curved to side, 9 1/2" (24.1 cm)




754370

TEALE UTERINE VULSELLUM

40 degree Curve Sideways
3 x 4 Teeth

- 754370 9 1/4" (23.4 cm)



182910

THOMS-ALLIS FORCEPS

- 182910 6 x 7 teeth, 8" (20.5 cm)



213250

BILLROTH TUMOR FORCEPS

- 213250 4 x 4 teeth, 10 1/4" (26 cm)



213140

BLAND, CERVICAL TRACTION

- 213140 9 1/8" (23.1 cm)

GYN Vulsellum *(continued)*



213150

HENROTIN VULSELLUM

213150 8 1/4" (21 cm)



213157

JACOB'S VULSELLUM

213155 Straight, 8 1/4" (21 cm)

213157 Curved sideways, 8 1/4" (21 cm)



213160

KELLY VULSELLUM

213160 Curved sideways, 9 1/4" (23.5 cm)



213175

M.G.H. VULSELLUM

213175 3 x 3 teeth, 6" (15.2 cm)



OB/GYN INSTRUMENTS



GYN Curettes



212905

HEANEY BIOPSY CURETTE

212905 9 3/4" (24.8 cm)



212910

KEVORKIAN-YOUNG BIOPSY CURETTE

212910 11 3/4" (29.8 cm)



212930

NOVAK BIOPSY CURETTE

212930 9 5/8" (24.4 cm)
diameter 1/4" (0.7 cm)



212925

RANDALL BIOPSY CURETTE

212925 9 5/8" (24.4 cm)



212753

SIMS UTERINE CURETTES

Sharp blades. Malleable shafts.
10 7/8" (27.6 cm) length

- 212749 Size 0
- 212750 Size 1
- 212751 Size 2
- 212752 Size 3
- 212753 Size 4
- 212754 Size 5
- 212755 Size 6




212784

THOMAS UTERINE CURETTES

Blunt blades. Malleable shafts.
10 7/8" (27.6 cm)

- 212780 Size 1
- 212781 Size 2
- 212782 Size 3
- 212783 Size 4
- 212784 Size 5
- 212785 Size 6

GYN Placenta and Sponge Forceps



121417

FOERSTER SPONGE FORCEPS

10 1/2" (26.7 cm)

- 121415 Straight, smooth jaws
- 121417 Straight, serrated jaws
- 121425 Curved, smooth jaws
- 121427 Curved, serrated jaws



212500

GYN Placenta and Sponge Forceps *(continued)*



212505

BOZEMAN UTERINE DRESSING FORCEPS

- 212500 Double-curved, one large ring. 10 1/4" (26 cm)
- 212505 Standard model 10 1/4" (26 cm)
- 212507 Straight model 10 1/4" (26 cm)



212507



770725

BARRETT PLACENTA FORCEPS 10 1/4" (26 cm)


- 770725 Small Jaw 28 x 14 mm



213050

DE LEE CERVICAL GRASPING FORCEPS

- 213050 10" (25.4 cm)



770731

IOWA MENGERT MEMBRANE FORCEPS 5 x 6 Teeth

- 770731 10" (25.4 cm)



INCHES

OB/GYN INSTRUMENTS



GYN Placenta and Sponge Forceps (continued)



213055

KELLY PLACENTA CLAMP

213055 12 1/2" (31.8 cm)




770730

POLK PLACENTA FORCEPS

9 3/4" (24.7 cm)
Fenestrated Jaw with Cross Serrations

770730 Straight




770405

WALLACE CAESAREAN FORCEPS

25 mm Jaw

770405 8" (20.3 cm)



GYN Biopsy Forceps



752528

WITTNER BIOPSY FORCEPS

Stem 9" (22.8 cm)
Toothed basket Jaw 8 x 5 mm

752526 Straight 11 1/2" (29.2 cm)

752528 Curved 11" (27.9 cm)



213036

KEVORKIAN-YOUNG CERVICAL BIOPSY FORCEPS

213036 10 5/8" (26.9 cm)




213011

GYN Biopsy Forceps (continued)

THOMAS-GAYLOR CERVICAL BIOPSY FORCEPS

- 213010 Straight, 9 1/2" (24.1 cm)
- 213011 Curved shanks, 9 1/2" (24.1 cm)



766522

SCHUBERT BIOPSY PUNCH

Curved Jaw 33 mm
Punch 5 x 8 mm

- 766520 10 3/4" (27.3 cm) Angled to Side
- 766522 11" (27.9 cm) Straight



213040

UNIVERSITY OF PENNSYLVANIA EQUINE UTERINE BIOPSY FORCEPS

- 213040 Biopsy Basket: 4 mm x 17 mm long
Shaft Working length: 24" (62 cm)
Overall length: 27.6" (70 cm)



214172

GYN Hysterectomy Forceps


HEANEY HYSTERECTOMY FORCEPS

Double serrated jaws with cross-groove
Length 8 7/16" (21.4 cm)

- 214170 Light jaws, single groove
- 214172 Heavy jaws, single groove
Length, 8 7/16" (21.4 cm)
- 214175 Heavy jaws, double groove.
Length, 8 7/16" (21.4 cm)



214175



OB/GYN INSTRUMENTS



214180

GYN Hysterectomy Forceps (continued)

HEANEY-BALLANTINE HYSTERECTOMY FORCEPS

Longitudinally serrated jaws. Single groove.
Length 8" (20.5 cm)

- 214180 Straight
- 214182 Curved



214185

LONG HYSTERECTOMY FORCEPS

- 214185 Short, powerful longitudinally serrated jaws with 1 x 2 distal teeth, 7 1/2" (19.2 cm)




214161

PHANEUF UTERINE ARTERY FORCEPS

Short heavy jaws, 1 x 2 distal teeth
Angled

- 214161 8 1/2" (21.5 cm)



182670

ROCHESTER-OCHSNER HYSTERECTOMY FORCEPS

1 x 2 distal teeth

- 182670 Curved, 8" (20.5 cm)



214230

WERTHEIM-CULLEN PEDICLE CLAMP

- 214230 Angular longitudinally serrated jaws. 8 5/8" (21.9 cm)

GYN Hysterectomy Forceps *(continued)*



762500

PHANEUF UTERINE ARTERY FORCEPS

8 1/4" (20.9 cm)
32 mm Jaw with Cross Serrations
1 x 2 Teeth 4 mm Tip


762500 Straight



638750

WECK HYSTERECTOMY FORCEPS

638750 Curved Jaw 10 1/4" (26 cm)
638752 Straight Jaw 10" (25.4 cm)



214215

DEBAKEY HYSTERECTOMY FORCEPS

214215 9 3/8" (23.8 cm)



213230

SOMER, UTERINE ELEVATING FORCEPS

213230 9" (22.8 cm)



214081

GYN Scissors

BRAUN EPISIOTOMY SCISSORS

Wide lower blade, angled to side.

214081 8 3/4" (22.2 cm)



INCHES

OB/GYN INSTRUMENTS



GYN Scissors (continued)



141512

KELLY FISTULA SCISSORS
Short, sharp pointed curved blades.

141512 6 1/4" (16 cm)



463463

MAYO UTERINE SCISSORS
8 3/4" (22.2 cm)

463463 Straight 53 mm blades
463464 Curved 54 mm blades



342242

NELSON UTERINE SCISSORS
Medium weight with rounded blades, blunt tips.

342242 10" (25.5 cm)



142235

TUNGSTEN CARBIDE SCISSORS

142235 9" (22.8 cm) Curved
142237 9" (22.8 cm) Straight



142237



511251

GYN Needle Holders

HEANEY NEEDLE HOLDER

Regular Jaw 22 x 3 mm

511252 8" (20.3 cm)



HEANEY NEEDLE HOLDER

Bulldog Jaw
Cross Serrations and Curved Jaw 15 x 3 mm

511251 8" (20.3 cm)



511250

HEANEY NEEDLE HOLDER

Tungsten Carbide Inserts
Curved Jaw with Cross Serrations 12 x 3 mm

511250 8" (20.3 cm)



510100

MAYO HEGAR NEEDLE HOLDER

Tungsten Carbide Inserts
Narrow Pattern

- 510100 6" (15.2 cm) Jaw 14 x 3 mm
- 510110 7" (17.7 cm) Jaw 14 x 3 mm
- 510120 8" (20.3 cm) Jaw 16 x 4 mm
- 510130 9" (22.8 cm) Jaw 20 x 2.5 mm
- 510140 11.75" (29.8 cm) Jaw 23 x 2.8 mm




508670

GYN Laser Instruments

AUVARD WEIGHTED VAGINAL SPECULUM

Ebonized for use with laser surgery.

- 508670 3" (7.5 cm) blade length
- 508671 3 3/8" (8.5 cm) blade length
- 508672 4" (10.1 cm) blade length



508604

BELLINA ANGLED SCISSORS

508604 9 1/2" (24 cm) working length



INCHES

OB/GYN INSTRUMENTS



GYN Laser Instruments (continued)




508603

BELLINA CERVICAL PUNCH SET

Set consists of two handles and three interchangeable punches.

- 508603 Punch Set with 5 1/8" (13 cm) and 15 1/2" (24 cm) handle lengths.




506483

BELLINA COAGULATING FORCEPS

Used to grasp and coagulate a bleeding vessel. The CO2 laser beam aimed at the ebonized forceps tip heats the tip, and the heat coagulates the vessel held in the jaws of the forceps.

- 506483 9 1/16" (23 cm) working length, serrated ebonized jaws



506482

BELLINA MICROSURGICAL BIOPSY FORCEPS

Double-action jaws

- 506482 9 1/16" (23 cm) working length



508600

BELLINA RULERS FOR LASER SURGERY

For determining depth and length of tissue for ablation purposes.

- 508600 Lateral ruler, graduated in mm up to 15 mm.
- 508601 Straight ruler, graduated in millimeters with a ridge indicating the 5 mm mark.



508601



508680

BOZEMAN GRASPING FORCEPS

Ebonized for use with laser surgery.

- 508680 Double curved, 10" (25.5 cm)
- 508681 Straight, 10" (25.5 cm)

GYN Laser Instruments *(continued)*

CERVICAL LASER HOOKS

Ebonized for use with laser surgery.

- 508605 Small, 2 mm hook, 6 1/2" (16.5 cm) working length
- 508606 Large, 4 mm hook, 7" (18 cm) working length

GELPI PERINEAL RETRACTOR

Ebonized for use with laser surgery.

- 508673 Single, sharp prongs with self-retaining ratchet mechanism.

GOODELL UTERINE DILATORS

Ebonized for use with laser surgery.

- 508685 Small, 10 5/8" (26.9 cm)
- 508686 Large, 13" (33 cm)

GYN LASER MIRRORS

With smoke evacuation channel

High-polish stainless steel mirror surface, 1" (25 mm) width.

- 508628 Angled handle Mirror, 6 1/2" (16.5 cm) working length
- 508630 Straight handle Mirror, 45 degree angled mirror, 7 3/8" (18.7 cm) working length

JACKSON VAGINAL RETRACTORS

Ebonized for use with laser surgery.

- 508674 Blade size 1 9/16" x 3" (4 cm x 7.5 cm)
- 508675 Blade size 1 9/16" x 3 1/2" (4 cm) x 9 cm)
- 508676 Blade size 1 9/16" x 4" (4 cm x 10 cm)

JARCHO TENACULUM

Ebonized for use with laser surgery.

- 508692 8" (20.5 cm)



508606



508673



508686




508630



508676



508692



OB/GYN INSTRUMENTS



508690

GYN Laser Instruments *(continued)*

KEVORKIAN-YOUNG CERVICAL BIOPSY PUNCH
Ebonized for use with laser surgery.

508690 9 3/4" (24.7 cm)



509210

LASER MIRRORS

Highly polished stainless steel front surface mirror
Suction to draw smoke and steam
Available in five sizes

509208 4 mm diameter mirror,
8 1/4" (21 cm) working length

509209 6 mm diameter mirror,
8 1/4" (21 cm) working length

509210 8 mm diameter mirror,
8 1/4" (21 cm) working length

509212 10 mm diameter mirror,
8 1/4" (21 cm) working length

509213 14 mm diameter mirror,
8 1/4" (21 cm) working length



509204

MICRO SUCTION TUBE

Ebonized for use with laser surgery.

509204 7 Fr., 8 1/2" (21.5 cm) working length



508602

MILAN-LONDON VACUUM TUBE

Ebonized for use with laser surgery.

Clips onto a standard metal Graves speculum for evacuation of smoke and plume.

508602



GYN Laser Instruments *(continued)*

PILLING LASER SPECULA

Inner blades are matte finished to help reduce reflection of the laser beam. External blades are polished stainless steel for ease of insertion. One smoke evacuation tube on the upper blade for rapid removal of steam and smoke from the ablated tissue. Available in three sizes.

- 508610 Small Laser Speculum
- 508612 Medium Laser Speculum
- 508614 Large Laser Speculum



508614

SCHROEDER TENACULUM

Ebonized for use with laser surgery.

- 508694 10" (25.4 cm)



508694

SCHUBERT CERVICAL BIOPSY PUNCH

Ebonized for use with laser surgery.

- 508691 Angled, 10 1/4" (26.2 cm) long, tapered bite



508691

THOMS-ALLIS VULSELLUM

Ebonized for use with laser surgery.

- 508696 6 x 7 teeth, 8" (20.5 cm)



508696

OB/GYN INSTRUMENTS

GYN Laser Instruments (continued)



429606

PILLING TITANIUM GYN LASER PROBE SET

Titanium dissipates heat of the laser beam

Universal Handle is interchangeable

Laser Probe Set includes:

429606 Sterilizing Tray, 10" (25.4 cm) x 10 1/2" (26.7 cm)

508640 Round-tapered probe, 6" (15.2 cm), 5 mm, 20 degree

508644 Round-tapered Probe, 6" (15.2 cm), 5 mm, 90 degree

508646 Flat Probe, 6" (15.2 cm), 5 mm, 20 degree

508650 Flat Probe, 6" (15.2 cm), 5 mm, 90 degree

508652 Universal Handle, 3" (7.6 cm)



508640



508644



508646



508650



508652

